

Humanized IgG Göttingen Minipigs for pre-clinical safety assessment of therapeutic antibodies

Genetically altered Göttingen Minipigs carrying a mini repertoire of human Ig- γ 1/ γ 4 heavy and the human κ light chain genes are available at Ellegaard Göttingen Minipigs A/S. The humanized Göttingen Minipigs show tolerance to a broad range of human antibodies, providing a novel model for safety testing and an important and viable alternative for toxicology studies of therapeutic antibodies in non-human primates.

Background

Therapeutic antibodies have recorded a huge increase in sales in the pharmaceutical market and have experienced an immense growth in recent years. The majority of FDA approved monoclonal antibodies (mAbs) are in the format of human or humanized IgG^{1,2}. mAbs are increasingly used for a variety of targets in cancer, diabetes, autoimmune, inflammatory, respiratory, ophthalmological, and infectious diseases.

Toxicological testing of antibodies in vivo can be challenging as most human therapeutic antibodies trigger xeno responses in preclinical animal models, resulting in rapid clearance of the drug or toxicities. This can be circumvented by using genetically altered animals that express the human protein and therefore recognize it as self. Humanized murine models have been generated to evaluate the in vivo stability of human mAbs³. However, data obtained in mice are not always translatable into humans in terms of application routes and pharmacokinetics. Furthermore, mice differ significantly from humans in general physiology, anatomy, and immune mechanisms, and the murine models are not ideal for long-term toxicological evaluations.

Non-human primates (NHPs) are the closest species to humans in biological terms, and NHPs are historically the primary choice for testing of biologics⁴. However, NHPs are not humanized and will, like other species, develop an anti-drug antibody after repeated dosing of human molecules.

In addition, the FDA strongly encourages the use of appropriate alternative models to NHPs when possible⁵.

Göttingen Minipigs have for decades been accepted as a non-rodent large animal species by national authorities including FDA and EMA⁶ and is the most commonly minipig breed in pharmaceutical research⁷. Göttingen Minipigs have many advantages for preclinical studies, including similarities in physiological and pathophysiological responses and high level of immunological similarity^{8,9}. Genetically altered Göttingen Minipigs are thus a viable alternative to NHPs in many areas of preclinical research^{10,11}.

Ideal model for safety testing

As demonstrated by Flisikowska et al*, the humanized Göttingen Minipigs' sensitivity and immunogenic response to therapeutic antibodies make them an ideal model for safety assessments of human or humanized IgG based molecules as well as prediction of possible side effects.

The amount of human IgG protein expressed is sufficient to induce and preserve immunological tolerance to a range of human IgG antibodies tested including bevacizumab, daratumumab, atezolizumab, and cergutuzumab*. Furthermore, the model has shown to fully reflect the difference in immunogenicity between these antibodies.

Humanized Göttingen Minipigs are bred in a fully AAALAC accredited barrier facility at the premises of Ellegaard Göttingen Minipigs A/S, Denmark, with an equally high health status as that of the standard Göttingen Minipigs. Trait inheritance has been confirmed. The humanized Göttingen Minipigs are fully immunocompetent and are commercially available in large uniform groups, sexually mature, and ready to use.

*Flisikowska T., Egli J., Flisikowski K. et al. A humanized minipig model for the toxicological testing of therapeutic recombinant antibodies. *Nat. Biomed. Eng* (2022). <https://doi.org/10.1038/s41551-022-00921-2>

For more information, please contact:

Ellegaard Göttingen Minipigs A/S
Soroe Landevej 302
4261 Dalmoose, Denmark
Tel.: +45 5818 5818
ellegaard@minipigs.dk



www.minipigs.dk



References

1. Lu RM, Hwang YC, Liu IJ, et al. Development of therapeutic antibodies for the treatment of diseases. *J Biomed Sci.* 2020;27(1):1.
2. Grilo AM, A. The increasingly human and profitable monoclonal antibody market. *Science and Society* 2019;37(1):9-16.
3. Bazin R, Boucher G, Monier G, Chevrier MV, S., Broly H, Lemieux R. Use of hu-IgG-SCID mice to evaluate the in vivo stability of human monoclonal IgG antibodies. *Journal of Immunological Methods* 1994;172(2):209-217.
4. Namdari R, Jones K, Chuang SS, et al. Species selection for nonclinical safety assessment of drug candidates: Examples of current industry practice. *Regul Toxicol Pharmacol.* 2021;126:105029.
5. US Food and Drug Administration. Nonclinical considerations for mitigating nonhuman primate supply constraints due to Covid 2022.
6. van der Laan JW, Brightwell J, McAnulty P, Ratky J, Stark C. Regulatory acceptability of the minipig in the development of pharmaceuticals, chemicals and other products. *J Pharmacol Toxicol Methods.* 2010;62(3):184-195.
7. Heining P, Ruyschaert T. The use of minipig in drug discovery and development: pros and cons of minipig selection and strategies to use as a preferred nonrodent species. *Toxicol Pathol.* 2016;44(3):467-473.
8. Pabst R. The pig as a model for immunology research. *Cell Tissue Res.* 2020;380(2):287-304.
9. Dawson H. Comparative assessment of the pig, mouse, and human genomes: A structural and functional analysis of genes involved in immunity. In: McAnulty, P.A., Dayan, A., Hastings, K.H., Ganderup, N.-C., editors. *The Minipig in Biomedical Research.* Boca Raton, FL: CRC Press. p. 321-341. 2011.
10. Flisikowska T, Kind A, Schnieke A. Genetically modified pigs to model human diseases. *J Appl Genet.* 2014;55(1):53-64.
11. Kalla D, Kind A, Schnieke A. Genetically Engineered Pigs to Study Cancer. *Int J Mol Sci.* 2020;21(2).